



Loops and Needles




- Meet semiquantitative standard for handling samples
- Surface treated for enhanced droplet adhesion
- Non-toxic (USP Class VI)
- Problem free inoculation with smooth loop surface
- Precision molded
- Calibration documents available for any batch upon request
- No flaming required



Nunc Needles

Polystyrene. Sterile



Cat. No.	253988	254399
Description	Yellow needle	Yellow needle
Units per pack/carton/case	50/1000/4000	12/600/2400



BIOGENERICA S.R.L.
 Via Pescara 7 Int. 23 - 95030 Mascali (CT) - Italia
 Tel. +39 095914678 - Fax +39 0952246205 - info@biogenerica.it
www.biogenerica.it

Nunc Loops

Polystyrene. Sterile



Cat. No.	251586	254437	253287	254410
Description	Blue loop	Blue loop	Clear loop	Clear loop
Total volume, µl	10	10	1	1
Units per pack/carton/case	50/1000/4000	12/600/2400	50/1000/4000	12/600/2400